

Rapid Screening and Confirmational Analysis of Residual Pesticides in Agricultural Samples by GC ECD/FPD and GC-MS/MS

Matt Lasater, Jim Edwards, Meredith Conoley, Jessie Butler

• Thermo Electron Corporation

Overview

Purpose: Develop and test an intelligent software to streamline the analysis of pesticides in a seed oil matrix.

Methods: An extract of seed oil was spiked with various levels of pesticides. The samples were analyzed on a TRACE[™] GC Ultra gas chromatograph (GC) configured with a tandem Electron Capture Detector/Flame Photometric Detector (ECD/FPD) for the analysis of chlorinated and organophosphorus pesticides. Any samples with a detected concentration of a target compound above a predetermined limit were automatically added to a list of samples to be injected on the Polaris Q mass spectrometer (MS). An ion trap MS was chosen because of its superior quantitation ability in matrix via the MS/MS experiment. A robotic TriPlus[™] liquid autosampler with a 150 sample tray was used to inject on either the inlet that led to the GC detectors or to the inlet that led into the MS, as directed by the Smart Screening software program.

Results: Of the 44 pesticides studied, 35 were sufficiently halogenated to be detected at a 1 pg on column load with the ECD. 15 were detectable at 10 pg column load using the FPD, and all of the pesticides were seen on the PolarisQ using a two stage MS experiment at a 1 pg level. The linear range of the data showed excellent correlation constants of greater than 0.99. The linear range was from 1 pg to 100 pg injected on column for the mass spectrometer and ECD and was 10 pg to 1000 pg for the FPD. 18 replicate injections were made at the 10 pg_{inlet} level on the PolarisQ. The median precision for the mass spectrometer was 8%.

Introduction

Pesticides generally fall into two groups. The first is a class of halogenated products, like p,p'-DDT which has 5 chlorine atoms per molecule, that have a strong response on the ECD. The second class is a set of organophosphorus compounds that are based on acetylcholinesterase inhibition. This class is not always halogenated, but it always has a phosphate moiety. These compounds are observed selectively using an FPD with a phosphorus specific photometric filter. The two GC detectors were operated in tandem with the ECD as the base detector, since it is a non-destructive detector. The various detectors used in this work are shown schematically in Figure 2. When the screening indicated the presence of a pesticide, the software then automatically added a confirming and quantitating run on the mass spectrometer. The second sequence was automatically launched with no user intervention.





FIGURE 1. Diagram of the electron capture detector (ECD), flame photometric detector (FPD), and the PolarisQ ion trap mass spectrometer (MS).



Methods:

The seed oil extracts were prepared by taking 0.2 g of sample, performing a solvent extraction, then cleaned by gel permeation chromatography (GPC). The resultant extract was reconstituted to 1 mL with iso-octane. This extract was then spiked with the pesticide mix at 1, 10, 100, and 1000 pg_{inlet}. All injections were 1.0 µL.

Mass Spectrometer

A diagram of what a multistage mass spectral experiment entails is shown in Figure 1. The first stage of the MS/MS experiment is the isolation of the precursor ion. The isolation is then followed by the second stage of MS, collision induced dissociation (CID), of that precursor ion to generate product ions, which are subsequently scanned out of the ion trap. The variables used in the isolation and fragmentation of the precursor ion, the trapping well depth parameter u_0 , and the CID voltage are found in Table 1. These parameters are compound specific and must be determined experimentally. Table 1 lists the precursor ions that were isolated and the dominant fragment ions that were used for quantitation in the MS/MS experiments. A comparison of the full scan, or single stage MS, and MS/MS experiment is found in Figure 3.

Gas Chromatograph

A GC method for the separation of 44 pesticides was developed using the tandem ECD/FPD to detect both chlorinated and organophosphorus pesticides in a single injection. This injection was made on a 100% dimethylsiloxane stationary phase column. The oven was held at 80°C for 1.0 minute followed by a ramp of 30°C/min to 150°C which was held for 10 minutes. A final ramp of 7°C/min reached the final temperature of 230°C, which was then held for 10 minutes. The injection was made in a 200°C splitless injector with a splitless time of 1.0 minute. The carrier gas for both columns was helium flowing at 1 cm/min. Next, a method was developed for the analysis of the same pesticides on the PolarisQ by MS/MS. The oven was held at 90°C for 1.5 minutes. The temperature was ramped at 30°C/min until 190°C was reached. This temperature was held for 10 minutes and then ramped at 9°C/min until 230°C was reached. This temperature was held for 10 minutes. The injection technique for the PolarisQ used programmable temperature vaporization. It was used in a splitless mode for 1.0 minute. The initial temperature of the injector was 80°C and was held for 6 seconds. The injector was then ramped at 10°C/sec to a final temperature of 250°C, where it was held for 1.5 min. The injector was cleaned between samples by holding the temperature at 250°C for 10 minutes. Table 1 lists the retention times of the pesticides for both of the chromatographic conditions.

FIGURE 3. Methoxychlor spectra in full scan mode and MS/MS mode.

